

2X Taq PreMix (Master Mix)

Description:

2X Taq Premix contains Parstous Taq DNA polymerase, reaction buffer, dNTPs mixture, protein stabilizer, and optimizes the convenience to use by adding sediment for electrophoresis and 2X solution of loading dye. In general, 2X Taq Premix shows no decline of activity compare with Taq DNA Polymerase, even in a room temperature. 2X Taq Premix is good for under 3 Kb of PCR products.

Kit storage:

This kit should be stored at -20°C. Unnecessary repeated freeze/thawing should be avoided. Under these conditions reagents are stable for one year from the date of production.

General Reaction Protocol:

1. Thaw 2X Taq Premix.
2. Prepare a master mix as following table.

Component	Volume	Final conc.
2X Taq Premix	10 µL	1X
Upstream Primer (10 pmol/ µL)	1 µL	0.5 pmoles/µL
Downstream Primer (10 pmol/ µL)	1 µL	0.5 pmoles/µL
Template DNA	Variable	10 fg~1 µg
Sterilized D.W.	Variable	-
Total Volume	20 µL	-

3. Mix the master mix and dispense appropriate volumes into PCR tubes. Centrifuge the reactions in a microcentrifuge for 10 seconds.
4. Perform PCR using your standard parameters (3-step cycling).

Cycle	Time	Temp °C
1	5 Min	94
25 ~35	30 Sec	94
	30 Sec	50~60
	30 ~60 Sec	72
1	5 Min	72

5. Separate the PCR products by agarose gel electrophoresis and visualize with EtBr or any other means.

* For PCR products longer than 3~4 Kb, use an extension time of approximately 1 min. per Kb DNA.

* A DNA fragment which is amplified by Taq DNA Polymerase has A overhang, and it enables you to do cloning by using T-vector.