

سیستم الکترونیک و نرم افزار رسم منحنی electronic system & Interface card

Descriptions:

Recording the magnetization and hysteresis curves of a ferromagnet. In a ferromagnet, the magnetic induction $B = \mu_r \cdot \mu_0 \cdot H$; $\mu_0 = 4\pi \cdot 10^{-7} \text{ Vs/Am}$: magnetic field constant reaches a saturation value B_s as the magnetic field H increases. The relative permeability μ_r of the ferromagnet depends on the magnetic field strength H , and also on the previous magnetic treatment of the ferromagnet. Thus, it is common to represent the magnetic induction B in the form of a hysteresis curve as function of rising and falling field strength H . The hysteresis curve differs from the magnetization curve, which begins at the origin of the coordinate system and can only be measured for completely demagnetized material. In this experiment, a current I_1 in the primary coil of a transformer which increases (or decreases) linearly over time generates the magnetic field strength $H = I_1 \cdot (N_1/L)$; L : efficient length of iron core, N_1 : number of windings of primary coil. The corresponding magnetic induction value B is obtained through integration of the the voltage U_2 induced in the secondary coil of a transformer: $B = (1/(A \cdot N_2)) \cdot U_2 \cdot dt$; A : crosssection of iron core, N_2 : Number of windings of secondary coil. The computer-assisted measurement system interface card is used to control the current and to record and evaluate the measured values. The aim of the experiment is determine the relative permeability μ_r in the magnetization curve and the hysteresis curve as a function of the magnetic field strength H

Specifications :

pc oscilloscop : External oscilloscop 20 MHz

signal generator : .1 HZ ...100KHZ

Warranty : One year

After sales service : Ten years

