

Kia One-Step qRT-PCR SuperMix

Cat. No. AQ211

Storage: at -20°C in dark for one year

Description

Kia One-Step qRT-PCR SuperMix is a one-step qRT-PCR kit which provides high sensitivity, high efficiency cDNA synthesis and qPCR amplification. RNA template and reverse gene specific primer (GSP) are used for first-strand cDNA synthesis, and then qPCR is performed with resulting cDNA and forward/reverse GSP. All the reactions including reverse transcription and qPCR are completed in one tube and in one reaction system.

Advantages

- High efficiency cDNA synthesis with One-Step Enzyme Mix and Green qPCR SuperMix, followed by PCR amplification with resulting cDNA. Simple procedure helps to minimize contamination
- High sensitivity, high specificity, accurate data.

Applications

- Multiple copy and low copy gene detection
- Viral RNA and trace RNA detection


Kit Contents:

Component	Amount
One-Step Enzyme Mix	40 μ l
Green qPCR SuperMix	1 ml
Passive Reference Dye (50x)	40 μ l
RNAse Free Water	1 ml


Reaction Components 20 μ l

Component	Volume	Final Concentration
RNA Template	1pg->1 μ g	as required
Forward GSP (10 μ M)	0.4 μ l	0.2 μ M
Reverse GSP (10 μ M)	0.4 μ l	0.2 μ M
Green qPCR SuperMix	10 μ l	1x
One-Step RT Enzyme Mix	0.4 μ l	-
Passive Reference Dye (50x) (optional)	0.4 μ l	1x
RNAse-free Water	Variable	-
Total volume	20 μ l	-

Thermal cycling program (two steps)

45°C 5 min
 94°C 30 sec
 94°C 5 sec
 60°C 30 sec  40-45 cycles
 Dissociation step

Thermal cycling program (three steps)

45°C 5 min
 94°C 30 sec
 94°C 5 sec
 50-60°C 15 sec
 72°C 10 sec  40-45 cycles
 Dissociation step

FOR RESEARCH USE ONLY

Two-step method is suitable for high specific qPCR.
Three-step method is suitable for high efficiency qPCR.

Notes

- High quality RNA template is recommended for use to ensure successful cDNA synthesis .
- This kit is only suitable for GSP, but unsuitable for first-strand cDNA synthesis using Oligo(dT) or Random Primer.

FOR RESEARCH USE ONLY