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Electromagnetic molten metal stirring device

Integration of a resistance furnace and an electromagnetic stirrer in this device can help researchers in metallurgy. Electromagnetic stirrers are in a group of non-contact stirrers. In this system, a magnetic coil is placed around the melt crucible and applying of three-phase altering voltage can cause magnetic field. It is clearly known that the magnetic field effects on an electrical conductor can induce the electric current in that. Similarly because of induction of magnetic flux - according to Fleming's right hand rule - the current is induced in the metal inside the crucible.

The current along with the electromagnetic field can develop the electromagnetic force (Lorentz) - according to Fleming's left hand rule - in the molten metal. Lorentz force is usually rotational and can move melt along the electromagnetic field direction. The coils can be designed for achieving to any desired stirring pattern.

Some features of the device:

- **Melt homogenization with no segregation**

In this method, the stirring is much stronger than the other methods and the melt stirs quickly and uniformly. Thus molten metal will be uniform and without any segregation.

- **Alloying possibility**

Due to the high ability of the system to homogenize the melt, it can be effective in alloying.

- **Solidification studies**

Stirring the metal during solidification is one of the interesting aspects of the solidification studies. In this process by stirring the melt during solidification, the dendrites are partially broken and the solidification front becomes unstable. It can consider different microstructures.

- **Possibility of continuous heating during stirring**

Researchers can determine the stirring time. Also maintaining the fluidity - which is achieved by the furnace - can cause better stirring.

- **Distribution of the solid phase in the melt**

Melting method is one of the methods of making matrix metal composites. In this method, reinforcement particles are added into the melts. In this process, stirring the melt plays an important role. Also melt should maintain its fluidity during stirring. Electromagnetic stirring system causes the separation of the second phase agglomerates and their uniform distribution by vortex forces.

- **Reduction of melting time and energy consumption**

Since the stirrer causes melt flow, heat transfer to the cold metal rises and it makes the charge to melt faster. Also due to the uniform heat distribution in

the molten metal that caused by stirring, increased local temperatures near the heat source are minimized. Thus the heat absorption from the heat source increases. Therefore it can decrease the melting time as well as energy consumption.

- **Oxidation reduction**

Shortening the melting time and no increase in local melting temperature can lower the oxidation.

- **Easy maintenance**

Since there is no contact between the devices and the melt in this system, maintaining the system becomes easier and just limited in system protection against crucible's heat and electronic protective arrangements.

- **Possibility of accurate temperature control**

If the melting temperature is not uniform, it would not be possible to measure the melt temperature accurately, even with the best measuring tools. Stirring process with electromagnetic non-contact technic results in a uniform temperature in all over the molten metal which can provide temperature reporting possibility.

- **Possibility of independent control of stirrer and furnace**

The furnace and stirrer parts can be separated into two independent parts which each part can work independently. Detaching the parts makes the cavity space larger for the stirrer.

- **Possibility of controlling the atmosphere**

The atmosphere inside crucible chamber can be protected via inert gases during melting and stirring. (this possibility is an extra option based on costumers demand)

Specifications:

Resistance furnace, capable for reaching 1200° C
 Inner diameter of the furnace /Stirrer for crucible location: 120 mm
 Allowed height of crucible in furnace/stirrer: 100 mm
 Inner diameter of the stirrer: 250mm
 Weight: 100 kg

The stirrer details:

| Parameter | Connection type (Y) | Connection type (Δ) |
|------------------------------|---------------------|---------------------|
| Number of phase | 3 | 3 |
| Voltage [V] | 400 | 400 |
| Frequency [Hz] | 50 | 50 |
| Current [A] | 30 | 70 |
| Power factor | 0.3 | 0.3 |
| Flux density (in center) [T] | 0.56 | 0.86 |
| Stirring Speed [rpm] | 3000 | 3000 |
| Insulation Class | F | F |