




Innovator In Spectroscopy Equipment

# DLA2031



## DELAY LINE AMPLIFIER MODEL DLA2031

www.cfp.co.ir



CATALOG



NIM MODULES



FREE DOWNLOAD

# DELAY LINE AMPLIFIER

## MODEL DLA2031

### Features

- Differential input
- Super fine gain control
- Excellent overload recovery
- Optimum timing capabilities
- Gaussian/triangular shaping
- Selectable integration time constants
- Ideal for neutron–gamma discrimination
- Excellent high counting rate performance
- Delay line shaping for energy and time spectroscopy
- High stability
- Splash resistant design
- Compact single width NIM
- Excellent temperature stability
- Wide operational temperature range
- Ultra reliable industry standard design
- Being robust and suitable for use in challenging conditions



OPTIMUM TIMING CAPABILITIES  
& EXCELLENT HIGH COUNTING RATE  
PERFORMANCE

## Description

The DLA2031 from Control Farayand Pasargad (CFP) is intended for energy and time spectroscopy with scintillation detectors. It can also be used with proportional counters, semi-conductor detectors, and position-sensitive proportional counters. Its delay line shaped output signal is particularly well suited for high counting rate and timing applications. This particular type of output signal offers a more rapid baseline recovery than is possible with semi gaussian shaping amplifiers. The model DLA2031 provides excellent timing capabilities, either for leading-edge or zero crossing timing techniques, particularly when it is used with a CFP model 2028 timing single channel analyzer. Double delay-line shaping exhibits less timing jitter when compared with either the classical RC-shaping network or active filter networks, primarily due to the fast rise time and fall time of the double-delay-line shaped output pulse.

## Specifications

### Input(s)

#### Input power

DLA2031 powered from a standard NIM bin and power supply

#### INPUT

BNC connector on front-panel accepts either positive or negative inputs with rise time of 10 to 1000ns and decay time of 25 to 2000 $\mu$ s;  $Z_{in}$  @1000 W, dc-coupled; linear maximum 3.3V; absolute maximum 20V.

### Output(s)

#### OUTPUT(UNIPOLAR)

Prompt or delayed with full-scale linear range from 0 to +10V; single-delay-line shaped; baseline restored level adjustable to  $\pm 1V$ ;  $Z_o < 1W$ , dc-coupled through front-panel BNC connector;  $Z_o = 93W$ , dc-coupled through rear-panel BNC connector.

Short-circuit protected.

#### OUTPUT(BIPOLAR)

Prompt output with positive lobe leading, double-delay-line shaped with full-scale linear range of 0 to 10V;  $Z_o < 1W$ , dc-coupled through front-panel BNC connector;  $Z_o = 93W$ , dc-coupled through rear-panel BNC connector.

Short-circuit protected.

#### PREAMP OUT

rear-panel standard CFP power connector for mating preamplifier.

### Control(s)

#### FINE GAIN

Front-panel single-turn potentiometer for continuously variable gain factor of X0.3 to X1

#### COARSE GAIN

Front-panel seven-position switch selects gain factors of X10, 20, 50, 100, 200, 500, and 1000.

#### INTEG

Front-panel slide switch selects an integration time constant of 0.04, 0.1, or 0.25 $\mu$ s. For 0.04 $\mu$ s setting amplifier rise time is  $< 75$  ns.

#### PZ ADJ

Front-panel potentiometer adjusts pole-zero cancellation for decay times from 25 $\mu$ s to  $\infty$ .

#### POS/NEG

Front-panel slide switch sets input circuit for either input polarity.

#### DC ADJ

Front-panel potentiometer adjusts the dc-level for single-delay-line shaped unipolar output pulses.

#### DELAY IN/OUT

Rear-panel slide switch selects either delayed (In) or prompt (Out) timing for unipolar output pulses. Delay is equal to the width of the unipolar output pulse.

### Indicator(s)

No indicator

### Performance

#### Gain range

7-position coarse gain selection from 10 through 1000 and single turn fine gain control from 0.3 through 1; total gain is the product of course and fine gain settings.

#### Shaping filter

Front-panel switch permits selection of integration time constant with  $\tau = 0.04, 0.1, \text{ or } 0.25\mu\text{s}$  (40, 100, or 250ns).

#### Integral nonlinearity

$\leq \pm 0.05\%$

#### Noise

$\leq 20\mu\text{V}$  rms referred to input using 0.25 $\mu$ s integrate and maximum gain of 1000;  $\leq 25\mu\text{V}$  for gain = 50;  $\leq 60\mu\text{V}$  for gain = 10

#### Crossover walk

For constant gain, walk  $< \pm 1\text{ns}$  for 20:1 dynamic range;  $< \pm 2\text{ns}$  for 50:1;  $< \pm 2.5\text{ns}$  for 100:1. Crossover shifts  $< \pm 4\text{ns}$  for any adjacent coarse gain switch settings.

#### Count rate stability

A pulser peak at 85% of analyzer range shifts  $< 0.2\%$  in the presence of 0 to 105 random counts/s from a Cs-137 source with its peak stored at 75% of the analyzer range.

#### Overload recovery

Bipolar recovers to within 2% of rated maximum output in  $< 5$  non-overloaded pulse widths from X500 overload; unipolar recovers in same time from X100 overload.

#### Delay line shaping

1 $\mu$ s. both delay lines have the same value.

**Temperature instability gain**

Gain:  $\pm 0.01\%/^{\circ}\text{C}$ , 0 to  $50^{\circ}\text{C}$


DC Level:  $\leq \pm 0.1\text{mV}/^{\circ}\text{C}$ , 0 to  $50^{\circ}\text{C}$ .

**Application**

- BIT (built-in test)
- Gamma-Gamma coincidence
- Pulse –height spectroscopy
- Also be used in simple pulse –counting
- Pulse –amplitude (energy) spectroscopy
- Scientific
- Radar receivers
- General laboratory usage
- Magnify the amplitude of the preamplifier output pulse

## Electrical and Mechanical

### Unit: DLA2031

Version type	NPHV2039-PRO	NPHV2039A-PRO	NPHV2039B-PRO	NPHV2039C-PRO	NPHV2039D-PRO	NPHV2039E-PRO
<b>Electrical</b>						
Power required	Its power from a NIM bin power supply. Required dc voltages and currents are +24V, 90mA; -24V, 90mA; +12V, 85mA and -12V, 75mA.					
<b>Physical</b>						
Dimensions (L x W x H)	220mm x 34mm x 248mm (L x W x H)					
Weight	1.216kg					
Mechanical	 <p style="text-align: right;">Unit: mm</p>					
<b>Environmental</b>						
Storage temperature	- 40°C to + 85°C					
Operating temperature	0°C to +50°C					
Relative humidity	< 80%					

## Software and user interface

The device doesn't have any software

## Ordering info

### DLA2031 Standard package includes

Part #	Image	Description
DLA2031 main		Delay line amplifier model DLA2031
ACCE2031001		CD user guide (1 Pack)
ACCE2031002		Box with foam insert
ACCE2031003*		Guaranty (one year)

\* =we stand behind our products. We guarantee your satisfaction in the quality of our instruments by providing a complete one-year warranty covering any defect of workmanship, material, and/or design. If our products do not perform, we will provide complete repair and/or replacement. for guaranty conditions, please refer to manual device (DLA2031- Manual)

### Optional accessories and services

Part #	Image	Description
ACCE2031004		Installation
ACCE2031005		Training
ACCE2031006**		Re-calibration (interval) services. 1year factory maintenance suggested, not required
ACCE2031011		BNC terminator 50Ω
ACCE2031012		RG58A/U, 50Ω cable with two BNC male plugs
ACCE2031013		Conn housing plug 50POS AMP connectors
ACCE2031014		Conn pin hood int 50pos panel MT
ACCE2031015		Guide pin 4-40
ACCE2031016		TE connectivity AMP connectors multimate, type II series pin
ACCE2031017		Bin guide pin

ACCE2031018



Guide socket

\*\* = The proper maintenance & calibration of your instruments is critical to ensure proper performance & accuracy. for Re-calibration (interval) services, please call with CFP company (021- 4604538)







Innovator in Spectroscopy Equipment



Unit 10, No 64, Vahedi  
(7) St, After Punak Sq,  
Ashrafi Esfahani Expy,  
Tehran, Iran



+98 (21) 46045383

w w w . c f p . c o . i r