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
7 Cleaning / Maintenance

All parts of the sampler can be cleaned with soap solution (then rinsed with distilled water and dried at room temperature) or disinfect with isopropanol .

All seals are maintenance-free .

Pistons should be lubricated lightly (with enclosed silicone grease) after cleaning or when new .

The sampler can be disassembled completely for cleaning and maintenance . See Fig .1 for individual parts .

- Hold pipette horizontally .
- Hold tip ejector button (3) with a thumb and pull out and pull off ejection sleeve (11) .
- Unscrew lower housing (10) of pipette by turning counterclockwise . If necessary , use hole A of pipette wrench , Fig . 4 (a) .
- Remove piston (5) and spring (6) .
- Unscrew clutch (5) with hole C of pipette wrench , Fig . 4 (b) .
- Take out nose cone (9) .
- Take out piston guide . Where appropriate , replace seal (8) .
- Remove residue in the nose cone :
Sampler 200-1000 Micro liter :(5)
Aspirate distilled water through the nose cone , e.g . With a water jet pump .
Dry at room temperature .

Sampler 5-100 Micro liters :
Push out reducing tube using the wire punch (accessory) , see Fig . 5 and replace .

- Swollen or damaged seals must be exchanged .

Assemble sampler by following above points in reverse order (Fig . 1) .

Note

- When inserting nose cone into lower housing , turn cone until it locks into position .
- Screw clutch into the lower housing of the sampler using the pipette wrench .
- Screw upper and lower housing together without using the wrench .
- Hold pipette horizontally . Push the ejection sleeve into the upper housing .

Caution :

After exchanging parts or other maintenance , always check that the pipette functions correctly . If a problem cannot be solved with the aid of the recommendations above , return the sampler to your local distributor or the manufacturer .

6 Troubleshooting table

<i>Error</i>	<i>Cause</i>	<i>Solution</i>
Droplets on the inside of the pipette tip .	Nonuniform wetting of the plastic .	- Attach new pipette tip .
Sampler drips , pipetted volume incorrect .	Tip is loose . Wrong pipette tip .	- Push tip on tightly . - Use original tip .
	Sampler leaks because :	
	- Piston contaminated	- Clean piston .
	- Seals damaged	-Exchange seals .
	- Nose cone loose .	-Tighten nose cone .
Control button jams , moves erratically .	Piston contaminated . Swelling of seal due to careless pipetting of organic solvents or caustic liquids .	- Clean piston . -Disassemble sampler and ventilate . Where appropriate , clean piston and seal with soap solutions , rinse with distilled water and lubricate lightly . Otherwise , change seal .

If the problem is still exist , return it to your local dealer or manufacturer for repair and calibration .

5 Testing the sampler

The sampler can be tested by weighing the volume of the bidistilled water pipetted on a sufficiently sensitive analytical balance .

Important :

The bidistilled water , weighing vessel , pipette and pipette tip must have reached the same temperature . Weighing should take place at 20 - 25 °c (constant to ± 0.5 °c) .

Avoid drafts !

To calculate the volume , divide the weight of the water by its density (at 20 °c : 0.9982) .

Sampler \leq 10 Micro liters

Testing is performed by aspirating the volume from a weighed vessel filled with water with a tip which has not been pre-wetted .

Sampler > 20 Micro liters

Bidistilled water is dispensed into a vessel from a pre-wetted tip . 5 measurements are sufficient for accuracy control and 15 measurements for precision control .

3.1 Special notes

Sampler > 10 Micro liter :

For utmost precision and accuracy , we recommend that liquid is aspirated and dispense two or three times when a new pipette tip is used :

- Aspirate and dispense liquid into tip several times .
- Empty the tip completely against the inside of the vessel outside the liquid (blow-out) .

Physical reasons :

- Pressure compensation of the system .
- Compensation of slight differences in temperature (Sec . 4) .
- Compensation of properties of the liquid . Pre-wetting of tip : wetting liquids (serum , detergent) from a thin film on the inside of the pipette tip . Therefore , the first volume dispensed would be too small .

When pipetting serum or solutions with a high viscosity , wait a few more seconds when aspirating and dispensing .

Sampler \leq 10 Micro liter :

The given values for accuracy and precision of small sample quantities can only be achieved when the sample is rinsed in another liquid .

The tip must not be pre-wetted during this process :

- Immerse tip containing sample into this liquid .
- Operate measuring stroke several times .
- Perform blow-out and slide tip up against inside of the vessel .
- Discard tip .

4 Operational limitations

- Differences in the temperature of the pipette and solution result in incorrect dispensing .
- In addition , volume errors may occur when pipetting liquids with a high vapor pressure .
- If there is no calibration facilities do not touch the calibration screw . (Fig3)

3 Mode of operation

The color of the control button indicates the color of the appropriate pipette tip .

If volumes > 10 Micro liter are pipetted , we recommend that liquid is aspirated and pipetted two or three times before dispensing when a new pipette tip is used . (See special note 3.1) .

Volumes ≤ 10 micro liter should be rinsed from an unwetted tip in another liquid (see special note 3.1)

Filling

- Attach appropriate pipette tip tightly (Pay attention to color code!)
- Press control button down to first stop (measuring stroke) , Fig . 2 (a) .
- Hold pipette vertically and immerse tip approx . 3 mm into the liquid .
- Let control button glide back slowly , Fig . 2 (b)
- Slide tip out of the liquid along the inside of the vessel .
- Wipe off any droplets with lint-free tissue . Ensure that no liquid is aspirated out of the tip .

Dispensing

- Hold tip at an angle against the inside of the vessel , Fig . 2 (c) .
- Press control button slowly down to first stop (measuring stroke) and wait approx . 1-3 seconds , depending on volume , Fig . 2 (c) .
- Press button down to second stop (blow-out) to empty tip completely , Fig . 2 (d) .
- Hold down control button . Slide tip up against the inside of the vessel .
- Let control button glide back .
- Eject the tip by pressing the separate tip ejector button , Fig . 2 (e) .

It is advisable to store the sampler in a pipette stand during breaks .

Never lay down the sampler with a filled tip !

1 Design principle

The sampler (piston stroke pipette) with separate tip ejection mechanism becomes a functional unit when an appropriate pipette tip is attached .

The liquid to be pipetted only becomes in contact with the pipette tip .

The material of the sampler is largely resistant to organic solvents. However , no liquid should penetrate into the sampler .

The buttons of the sampler have the following functions:

control button :


1st stop (measuring stroke) : the volume is aspirated or dispensed .
2nd stop (blow - out) : the remaining liquid in the tip is dispensed .

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

micropipetting system Instruction manual

e-mail: info@labtronco.com
www.labtronco.com

labtron