Taligene Pars

# Paraoxonase 1 Activity Assay Kit

Catalog No: TGP-PON-2011  
96T

Taligene Pars Company No. 16,  
Commercialization Complex, Isfahan Science  
and Technology Town, Isfahan, Iran.

Taligene@yahoo.com  
Phone +98 3133932373  
+98 9130390282  
Fax +983133932374

<http://www.taligene.com>

## I. Introduction

Paraoxonase 1 (PON1, EC 3.1.8.1) is an HDL-associated esterase/lactonase enzyme mainly synthesized in the liver and released into the bloodstream. PON1 hydrolysis highly toxic organophosphorus compounds such as the pesticide, paraoxon and various nerve agents used as chemical weapons. Additionally PON1 has protective effect against oxidative damage has anti-atherogenic properties and protects serum HDL and LDL particles against lipid peroxidation. PON1 inhibits N-homocysteinylolation of LDL-associated proteins by hydrolyzing the highly reactive pro-oxidant homocysteine thiolactone. Measurement of serum PON1 activity has been proposed as an antioxidant capacity index and biomarker of liver and cardiovascular disease.

## II. Principle of Assay

Taligene's Paraoxonase 1 Activity Assay Kit enables rapid measurement of PON1 activity by only one assay reagent. In this assay, PON1 catalyzes the cleavage of paraoxon resulting in p-nitrophenol formation. The rate of formation of p-nitrophenol is measured by monitoring the increase in absorbance at 412 nm at 25°C.


Fig. 1. Reaction of PON 1 (1).

## III. Kit Components

The working reagent consists of Tris/HCl buffer, pH 8.0, containing  $CaCl_2$ , paraoxon as the substrate and stabilizer. Reagent is ready to use.

## IV. Storage Conditions

Store components at 2-8 °C until their expiration dates.

## V. Samples

Blood serum, heparinized plasma, seminal plasma, cell lysates, and tissue homogenates can be used as sample.

## VI. Pre-assay preparation

- Sample preparation
- **Serum:** Use a serum separator tube (SST) and

allow samples to clot for two hours at room temperature or overnight at 4°C before centrifugation for 15 minutes at 1000 ×g. Remove serum and assay immediately or aliquot and store samples at -20°C or -80°C. Avoid repeated freeze-thaw cycles.

–**Plasma:** Collect plasma using heparin as an anticoagulant. Centrifuge for 15 minutes at 1000 x g. Assay immediately or aliquot and store samples at -20°C or -80°C. Avoid repeated freeze-thaw cycles.

–**Cell Culture Supernatant:** Remove particulates by centrifugation for 15 minutes at 1000 x g, 2 - 8°C and assay immediately or aliquot and store samples at -20°C or -80°C. Avoid repeated freeze-thaw cycles.

–**Tissue Homogenates:** You'd better get detailed references from other literatures before assay aiming at different tissue types. In general, tissue pieces should be weighed and then homogenized in PBS (the volume depends on the weight of the tissue) with a glass homogenizer on ice. The homogenates are then centrifuged for 5 minutes at 5000×g to get the supernatant.

**Note:** As PON1 is a strongly  $Ca^{2+}$  dependent enzyme, heparinized plasma samples should be used. Plasma specimens collected with EDTA or other  $Ca^{2+}$  chelating anticoagulants may exhibit reduced PON1 activity.

## VII. Assay procedure

**Step 1:** Turn on spectrophotometer and allow instrument to initialize for 15 minutes.

**Step 2:** Clean cuvettes with alcohol. Rinse well with deionized water.

**Step 3:** Zero spectrophotometer at 412 nm with 1000 µl working reagent (Blank).

**Step 4:** Pipette the following into the cuvette:

- 1000 µl working reagent
- 25 µl sample

**Note:** The reaction should be perfectly linear and there should be no plateauing. These would indicate the need to dilute or concentration samples further. In manual working, the volumes of the sample could be 5-100µl.

**Step 5:** Pipette up and down to mix thoroughly, while avoiding bubble formation.

**Step 6:** Place cuvette in the correct position in the

spectrophotometer. Record the change in  $A_{412}$  for 2 min.

## VIII. Calculations

One unit of paraoxonase activity is the amount of enzyme that generates 1 µmole of p-nitrophenol per min at 25°C and pH 8.

Molar extinction coefficient for p-nitrophenol is  $18290 M^{-1} cm^{-1} = 0.018290 \mu M^{-1} cm^{-1}$

A theoretical unique factor is determined to convert change in absorbance per minute ( $\Delta A/min$ ) to the corresponding units of enzyme activity. This factor is calculated using the following equation:

$U/L = \Delta A/minute \times F$ ; where F= factor

$F = (TV/SV) / 0.018290$  where;

TV = Total Volume in µl

SV = Sample Volume in µl

0.018290 = micro molar extinction coefficient

This factor can be programmed into the spectrophotometer and the machine directly converts the change in absorbance at 412 nm ( $\Delta A/min$ ) to activity in U/L.

## IX. Example of Results

PON1 activity was assayed at 25°C with a cuvette path length of 1 cm using the assay procedure above. Change in absorbance was recorded every 20 seconds up to 2 minutes.

$U/L = \Delta A/minute \times F$

initial absorbance=0.0

final absorbance = 0.6

$\Delta A/min = 0.3$

Paraoxonase activity  $U/L =$

$(0.3) \times (1025/25) / 0.01829 = 672.4$

Note: Results for diluted samples must be multiplied by their dilution factor.

**Warning:** Wear gloves and masks while performing the task.

## X. References

1. Mackness B, Mackness M, Aviram M., Paragh G. The Paraoxonases: Their Role in Disease Development and Xenobiotic Metabolism. Dordrecht; Netherland, Springer, 2008:3-260.
2. Bełtowski J, Wójcicka G, Jamroz A. Leptin decreases plasma paraoxonase 1 (PON1) activity and induces oxidative stress: the possible novel mechanism for proatherogenic effect of chronic hyperleptinemia. Atherosclerosis 2003;170(1):21-9.

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۱. Mackness B, Mackness M, Aviram M., Paragh G. The Paraoxonases: Their Role in Disease Development and Xenobiotic Metabolism. Dordrecht; Netherland, Springer, 2008;3-260.

۲. Beltowski J, Wójcicka G, Jamroz A. Leptin decreases plasma paraoxonase 1 (PON1) activity and induces oxidative stress: the possible novel mechanism for proatherogenic effect of chronic hyperleptinemia. Atherosclerosis 2003;170(1):21-9.