


PRM2415A

SIMPLE










PORTABLE RADIATION MONITOR MODEL PRM2415A SIMPLE













DOSIMETERS

FREE DOWNLOAD

CATALOG

PORTABLE RADIATION MONITOR MODEL PRM2415A SIMPLE


Features

- Detects all Gamma & X-Ray radiation
- Measurement Count:(0.01cpm~300Kcpm)
- Dose Range: (0.1uSv/h~10mSv/h (@Cs137))
- Compact design with Low weight
- Cost effective Semiconductor Sensor
- 2 versions STD type (120keV-3MeV) and PRO type (60keV-3MeV) gamma energy
- Light and Sound alarm (PRO version only)
- Detect the presence of a source of radiation During air travel or during medical procedures (X-Rays, CT scans, etc.)

- Monitor counts/interval (rate)
- High sensitivity and wide measurement range
- Monitor variation in background radiation at different elevations.
- Monitor radioactivity in the environment over long periods of time.
- Monitor counts per interval (rate) from a hard beta or gamma radiation source as a function of the distance between the source and the Radiation Monitor

X- and Gamma-Rays





- 2 Thin plates made of wood, aluminum, etc
- (3) Lead, iron, and other thick metal plates

1 Paper

Description

The PRM2415A from CFP is a highly sensitive radiation sensor designed for the embedded systems market. Capable of detecting Gamma and Beta radiation, it has a simple pulsed output that can be used with any interface. The PRM2415A has an onboard DC boost circuit, so the board can be supplied with a friendly 5V. Using only 31mW (6mA @ 5V), it is very low power. CFP has a handful of documents and example sample to get you up and running. They have also written a Windows and Android example program written in C# to output graphs to a computer using a USB as the reader. Additionally, the counter comes with an optional enclosure and 3.5mm cable that allow you to connect the counter to some mobile phones. Portable combined multipurpose instrument, which can be equipped with different external detection units for various applications as applicable.

The Radiation Monitor senses ionizing radiation by means of a SEM. The tube is fully enclosed inside the instrument. When ionizing radiation or a particle strikes the SEM surface, it is sensed electronically and monitored by a computer, or by a flashing count light. Radiation is measured in counts in a time interval,

as configured in data-collection software. About 5 to 25 counts at random intervals (depending on location and altitude) can be expected every minute from naturally occurring background radiation.

Further Tips for Monitoring Radiation

To measure gamma and X-rays, hold the tip of the Radiation Monitor toward the source of radiation. Low-energy gamma radiation (10-40KeV) cannot penetrate the side of the PRM, but may be detected through the window. To detect hard beta radiation, point the window toward the source of radiation. Beta radiation has a longer range through air, but can usually be shielded (e.g., by a few millimeters of aluminum). High energy beta particles may be monitored through the back of the case. Some isotopes detected relatively well are cesium-137, cobalt-60, technetium-99m, phosphorus-32, and strontium-90. Some types of radiation are very difficult or impossible for this PRM to detect. Beta emissions from tritium are too weak to detect using the Radiation Monitor. Gamma emissions from Americium-241 are too weak to detect using the Radiation Monitor well.

Specifications					
	UNIT:2415A				
ТҮРЕ	SIMPLE				
Version	STD	PRO			
	Inputs/	Outputs			
Port	USB	USB			
	Controls/Co	mmunication			
USB	Set and read all data	Set and read all data			
	Perfo	rmance			
Energy range	120keV-3MeV	60keV-3MeV			
Battery type	Not included	Not included			
Dose range	0.1-10000µSv/h	0.1-10000µSv/h			
Accuracy	+/- 20% (Cs-137)	+/- 10% (Cs-137)			
Power	0.1W	0.1W			
Protection class	IP40	IP40			
	Continuou	us run time			
In normal conditions	≥1000 h	≥1000 h			
In economy mode	≥5000 h ≥5000 h				
Buzzer	Not included	Available			
		cators			
Ring	Green color (Safe) Red Color (Alarm)	Green color (Safe) Red Color (Alarm)			

Application

- Research activities
- Nuclear medicine and industry
- Lantern Mantles
- Customs and border control
- Civil protection, security and police
- Emergency rescue service
- Individual radiation and protection
- Radioecology and Radiology
- Belding protection and emergency service
- Medical and health physics
- Accelerating installation
- Search X-ray apparatus
- Background Radiation
- Detects all Gamma & X-Ray radiation: Check radiation levels during air travel or during
- Counts/Interval vs. Distance and Shielding Studies
- Half-Life Determination (Counts/Interval vs. Time)
- medical procedures (X-Rays, CT scans, etc.)
- Dose monitoring of population
- X-Ray crystallography and X-Ray fluorescence spectroscopy, electronic microscopy

Electrical and Mechanical

VERSION

STD

- PRO
- Power required
- 0.1W
- weight
- 98 gr

Dimension

31mm x 137mm

Storage temperature

10°C...+40°C

Operating temperature

0°C...+50°C




STD Version



Software

This sensor can be used with one of the below interfaces and the following data-collection software.



Software main window

Print Preview				
Save As PDF	Print			
	Date: 1396-04-07 15:38:48 Device: PRM2415A SN: 181AITW9	Analysis report	کنتل فارند باسترگاد. کنتل فارند باسترگاد	
	Info: Analysis ID: 000003 Analysis mode: Gamma PreSet Timer: 1 Min	Dose Ratio Constant (Cs-137): 1CPM = 0.0017 Dose Ratio Constant (Co-60): 1CPM = 0.0019 Dose Ratio Constant (Sr-90): 1CPM = 0.0036	9mR/h = 0.0199µSv/h	
	Result: Y: 0.23 CPS 0.25 µ5∨/h 0.02 mR/h			
	Spectrum:			
	o Description:	20 40 Time(Sec) Gamma Spectrum	60	
	Responsible Expert:	Signature:		
Page 1 of 1				🔄 🗐 🗎 73% 🖨

Ordering info

tandard package includes PRO Version				
Part #	Image	Description		
PRM2415ASIMPLE-PRO	8	RRM 2415A-PRO mountable version		
ACCE2415A001-Simple		Phone jack 3.5mm cable		
ACCE2415A002-PRO		CD User guide (1 Pack)		
ACCE2415A003-PRO	GEZ	box with foam inserted (1 Pack)		
Standard package includes STD Version				
Part #	Image	Description		
PRM2415ASIMPLE-STD	8	RRM 2415A-STD mountable version		
ACCE2415A001-Simple		Phone jack 3.5mm cable		
ACCE2415A002-STD		CD User guide (1 Pack)		
ACCE2415A003-STD	GCP	box with foam inserted (1 Pack)		
Optional accessories and services				
Part #	Image	Description		
ACCE2415A003		installation		
ACCE2415A004		Training		



Innovator in Spectroscopy Equipment



Unit 10, No 64,Vahedi (7)St, After Punak Sq, Ashrafi Esfahani Expy, Tehran, Iran

+98 (21) 46045383

www.cfp.co.ir