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# mbcR 190 RQ (V3.0)

---

**mbcr 190 RQ** kit is intended for the detection of BCR-ABL (p190) transcripts in peripheral blood and BCR-ABL% calculation in patients undergoing therapy. This kit is designed for use with Rotor Gene 6000/Q or StepOne/StepOnePlus instrument. This kit is for research use only!

***Important Note:*** *This kit doesn't provide reagents for RNA extraction or cDNA synthesis!*

## 1. Introduction

Philadelphia chromosome is an abnormality resulted from 9;22 translocation. Consequently ABL proto-oncogene on chromosome 9 is fused with BCR gene on chromosome 22. This fusion produces BCR-ABL protein, mostly 210kDa (b2a2 or b3a2) or 190kDa (e1a2), with constitutively active tyrosine kinase activity promoting cell proliferation and inhibition of apoptosis. The fusion gene transcript is detectable in about 95% of CML patients and some cases of ALL. Also, serial monitoring of patients for identifying and measuring BCR-ABL transcripts provides more precise assessment of response to specific therapies and prediction of those in higher risk of disease progression.

This kit provides reagents for both detection of BCR-ABL transcripts (e1-a2 or p190) as well as BCR-ABL% calculation after therapy.

Currently Real-Time PCR provides the highest sensitivity and widest dynamic range among other methods for detection of BCR-ABL transcripts. In this method application of fluorescent dye labeled probes allows detection of amplified

# mcr 190 RQ (V3.0)

---

product. Analysis of fluorescent kinetics also leads to detection of the target sequence in the reaction without requiring post-amplification analysis, reducing the possibility of contamination with the PCR product.

mcr 190 RQ kit provides a ready-to-use Real-Time PCR system for detection of BCR-ABL (p190) transcripts with RotorGene 6000/Q (Qiagen) or StepOne/StepOnePlus (Applied Biosystems).

## 2. Kit Contents

The kit contains a manual, a CD with RotorGene and StepOne templates and following reagents:

Label	Content	Quantity
<b>mcr Mix*</b>	Master mix for BCR-ABL	500 µl
<b>ABL Mix*</b>	Master mix for ABL	500 µl
<b>m1</b>	mcr Standard 1: 100,000 copy/µl	100 µl
<b>m2</b>	mcr Standard 2: 10,000 copy/µl	100 µl
<b>m3</b>	mcr Standard 3: 1,000 copy/µl	100 µl
<b>m4</b>	mcr Standard 4: 100 copy/µl	100 µl
<b>m5</b>	mcr Standard 5: 10 copy/µl	100 µl
<b>A1</b>	ABL Standard 1: 100,000 copy/µl	100 µl
<b>A2</b>	ABL Standard 2: 10,000 copy/µl	100 µl
<b>A3</b>	ABL Standard 3: 1,000 copy/µl	100 µl
<b>A4</b>	ABL Standard 3: 100 copy/µl	100 µl
<b>Water</b>	PCR Grade Water	150 µl

\* 1, 2 and 4 tubes for 24, 48 and 96 reaction kits respectively.

### 3. Storage and Stability

The kit components should be shipped and stored at  $-20^{\circ}\text{C}$  and are stable until the expiry date mentioned. Avoid repeated freeze-thaws especially for PCR MIX more than few times to prevent reduced sensitivity.

### 4. General Precautions

In order to prevent false results, always pay attention to the following points:

- **Treat all samples as potentially infectious.**
- Within the pre-PCR work area assign three separate spaces for: a) Sample storage and extraction, b) Reagent preparation where the master-mix is aliquoted into tubes and c) Reaction preparation area for addition of templates to the tubes.
- Always wipe the working surfaces with 70% Ethanol before and after work.
- Thaw kit components completely, mix by flickering followed by a quick spin and store on crushed ice after.
- Keep PCR Mix tube at  $-20^{\circ}\text{C}$  at all times. Take it out just before use and return it to freezer immediately after.
- Do not place 0.2ml PCR tubes on crushed ice. Use cooling blocks instead.

### 5. Additionally Required Materials

- Real-Time PCR machine and accessory computer
- Refrigerated table top microtube centrifuge (for RNA extraction)

- Vortex mixer
- Adjustable pipetters
- RNA extraction kit
- cDNA synthesis kit
- Nuclease free filtered tips
- Tubes or strips and strip caps
- Disposable powder-free gloves

## **6. Specimen, storage and transport**

Peripheral blood should be collected in sterile condition in proper and sterile tubes. We recommend EDTA or Citrate as anticoagulant. Whole blood should be shipped and stored at +4°C (stable for 72 hrs).

RNA can directly be extracted from whole blood. Alternately to increase sensitivity buffy coat can be used. For optimum results a sample should include 5 to 10 million WBC per 150µl.

Whole blood or buffy coat can be stored at +4°C for three days. Otherwise should be aliquoted and stored at -70°C which is stable for few months.

## **7. RNA isolation**

RNA isolation can be performed using different kits from various manufacturers. We recommend using TriPure isolation Reagent (Cat. no. 1667157, Roche Applied Science, Mannheim, Germany). We also recommend avoiding silica gel or column based kits for total RNA extraction.

## 8. cDNA synthesis

1µg of total RNA is required and should be reverse transcribed to cDNA using random hexamers. Different kits are available in the market for this purpose.

Dilute prepared cDNA 2.5x with nuclease free water. For example to 20µl of cDNA add 30µl of nuclease free water.

## 9. PCR Protocol

Thaw the reagents on ice completely followed by brief mixing and a quick spin.

Each sample should be examined for both BCR-ABL fusion (p190) gene and for ABL control gene in two separate set of reactions. In BCR-ABL set consider 1 tube for each sample as well as 6 tubes for the 5 standards (m1 to m5) and “no template control” (NTC). In ABL set consider 1 tube for each sample and 5 tubes for the 4 standards (A1 to A4) and NTC. Place required number of tubes on pre-cooled (0-4°C) cooling block or loading block.

**Pipette 20µl of mbc190 Mix to the first series of tubes and 20µl of ABL Mix to each tube of the second group. Continue by adding 5µl of cDNA or standard or control to each tube.**

Cap the tubes/Strips and visually inspect to make sure all are capped securely. Place tubes in the machine. If using StepOne/StepOnePlus instrument, spin tubes/strips briefly before loading tubes on the block. If using RotorGene attach the locking ring.

## **10. Programming RotorGene**

*Before you start the machine make sure you have attached the locking ring on the rotor!*

Open the CD provided in the kit and double click on M-BCR 0.2 (if 0.2ml tubes are used) or on M-BCR Strip (if 0.1ml tubes are used) to open the program. Click on Start button (Green button on the top menu). On the pop up window click start again and save the run file.

Edit sample names on both mbcR and ABL pages. Remember that tubes containing mbcR mix should only be named in mbcR page and tubes containing ABL Mix should only be named in ABL page.

Make sure in the "Type" column, all the standards have been defined as "standard" and the relative concentrations have been entered. Patient samples should be defined as "unknown" and no template control as "NTC" respectively.

## **11. Programming StepOne**

Open the StepOne software (V 2.\*). On the Set Up menu click on Template and select the file on CD provided with the kit. Click on Plate Setup. One negative control, 5 standards for mbcR, 4 standards for ABL and few samples are defined. You may change plate set up using right click options (copy, past, clear). You may also add or remove samples on "Define Targets and Samples" menu. When finished, click on Start Run and save the experiment on desired location. Instrument will start shortly.

## 12. Programming Other Machines

If you apply this kit to other Real-Time PCR machines, program it according to the following table:

Step	Temperature and time	Cycles
1	95C x 10 min	1
2	95C x 15 sec	45
	60C x 60 sec	

Fluorescence should be collected at 60C for FAM and VIC/HEX dyes.

Both of mbr Mix and ABL Mix contain ROX.

## 13. Data Analysis: RotorGene

Before analyzing results, make sure in the sample menu all the standards have been defined as "standard" and the relative concentrations have been entered. Patient samples should be defined as "unknown" and Negative control or no template control as "Negative Control" or "NTC" respectively.

Analyze data according to manufacturer recommendations. Perform quantitative analysis for both **mbr (Green channel)** and **ABL (Yellow channel)**. Briefly, click on analysis menu and then under Quantitation tab double click on cycling A. Green. Close the pop up window and manually set threshold at 0.05. Repeat above for Yellow channel and set the threshold on 0.05.

Figures 1 and 2 represent typical graphs for RotorGene machine.

# mbcR 190 RQ (V3.0)


Figure 1. Typical mbcR Graph in Green Channel for RotorGene


Figure 2. Typical ABL Graph in Yellow Channel for RotorGene

Consider following points when analyzing:

**Note that calculated CT is reliable only when there is a sigmoid graph with clear logarithmic phase. In the absence of sigmoid graph, CT must be ignored to prevent false results.**

- A sample is **Positive** if generates a sigmoid graph in mbcR (Green) channel with CT of 20-40 and also a

# mbc190 RQ (V3.0)

---

sigmoid graph in ABL (Yellow) channel with CT of 20-30.

- A sample is **Negative** if does not generate sigmoid graph for mbc190 (Green) channel while generates a sigmoid graph for ABL (Yellow) channel with CT of 20-27.
- Results are **Inconclusive** and the test should be repeated if a sample generates a sigmoid graph in ABL (Yellow) channel with a CT higher than 27 but fails to produce a sigmoid graph in mbc190/Green channel.

*Note: All patient samples should generate a sigmoid graph with a CT of 27 or less in ABL/Yellow channel (with threshold of 0.05). A CT greater than 27 usually happens if not enough cells have been extracted or less than 100ug RNA has been used. ABL/VIC CT higher than 27 reduces the sensitivity of the test.*

## 14. Data Analysis: StepOne

Analyze data according to manufacturer recommendations. Briefly, click on Analyze and set the threshold for both **mbc190/FAM** and **ABL/VIC** on 0.1.

Figures 3 and 4 represent typical graphs for StepOne machine.


Figure 3. Typical mbr Graph in FAM Channel in StepOne


Fig 4. Typical ABL Graph in VIC Channel in StepOne.

Consider following points when analyzing:

**Note that calculated CT is reliable only when there is a sigmoid graph with clear logarithmic phase. In the absence of sigmoid graph, CT must be ignored to prevent false results.**

- A sample is **Positive** if generates a sigmoid graph in mbr/FAM channel with CT of 20-40 and also a sigmoid graph in ABL/VIC channel with CT of 20-30.
- A sample is **Negative** if does not generate a sigmoid graph for mbr/FAM channel while generates a sigmoid graph for ABL/VIC channel with CT of 20-27.
- Results are **Inconclusive** and the test should be repeated if a sample generates a sigmoid graph in ABL/VIC channel with a CT higher than 27 but fails to

## mbcR 190 RQ (V3.0)

---

produce a sigmoid graph in mbcR/FAM channel. ABL/VIC CT higher than 27 reduce the sensitivity of the test.

*Note: All patient samples should generate a CT of 27 or less in ABL/VIC channel (with threshold of 0.05). A CT greater than 27 usually happens if not enough cells have been extracted or less than 100ug RNA has been used.*

### **15. BCR-ABL% calculation**

To assess the response to therapy, BCR-ABL% value for each patient can be calculated. This kit uses NCN method for this purpose (*Beillard E. 2003, Leukemia 17:2474*). In this method BCR-ABL% value is the BCR-ABL expression (titer) normalized by the ABL expression (titer) and then multiplied by 100%.