



DNase I

Cat No: DN9061 Volume: 1ml

Contents: DNase I (5U/ μ l) 1ml
10x Reaction buffer with mgcl₂ 1ml

Store at: -20 °C

Description

DNase I, (RNase free) is an endonuclease that nonspecifically cleaves DNA to release di-, tri-, and oligonucleotide products with 5'- phosphorylated and 3'- hydroxylated ends. DNase I acts on single- and double-stranded DNA, chromatin and RNA:DNA hybrids.

Features

Recombinant enzyme.

Purified for non-animal host with a lower level of intrinsic RNases.

Applications

Degradation of DNA template in transcription reactions

Removal of contaminating Genomic DNA From RNA sample

DNase I footprinting

Nick translation

Definition of activity unit

One unit of the enzyme completely degrades 1 μ g of plasmid DNA in 10 min at 37 °C. Enzyme activity is assayed in following mixture: 40 mM (PH 8.0), 10 mM MgSO₄, 1mM CaCl₂, 1 μ g of pUC19 DNA.

Source

AE Coli strain that carries an MBP fusion clone of Bovine Pancreatic DNase I.

Molecular Weight

29 KDa monomer

Quality control

دفتر فروش: تهران، خیابان انقلاب اسلامی، خیابان دوازده فروردین جنوبی، انتیتو پاستور ایران، مرکز رشد زیست فناوری

تلفن: 021 2508 2505, 021 5543 3092 , 021 6411

واحد تولید و تحقیقات: تهران، خیابان انقلاب اسلامی، خیابان دوازده فروردین جنوبی، انتیتو پاستور ایران، مرکز رشد زیست فناوری

Email:v.gene.pars@gmail.com, info@vgenepars.ir Web: www.vgenepars.ir



The absence of ribonucleases confirmed by appropriate quality test. Functionally tested for digestion of template DNA after in vitro transcription.

Storage buffer

The enzyme is supplied in: 10 mM Tris-HCL (PH 7.5), 10 mM CaCl₂, 10 mM MgCl₂, 50% (v/v) glycerol.

10xReaction Buffer with MgCl₂

100 mM Tris-HCL (PH 7.5 at 25°C), 25 mM MgCl₂, 1 mM CaCl₂.

Inhibition and Inactivation

Inhibitors: metal chelators, transition metal (e.g., Zn) in millimolar concentrations, SDS (even at concentrations less than 0.1%), reducing agent (DTT and meta- mercaptoethanol), ionic strength above 50-100 mM.

Inactivated by heating at 65°C for 10 min in the presence of EGTA or EDTA (use at 1 mol of EGTA/EDTA per mol of Mn²⁺/Mg²⁺).

Note

DNase I is sensitive to physical denaturation. Mix gently by inverting the tube. Do not vortex.