



# **High Pure PCR & Gel Cleanup Kit**

**Catalogue Number: N-0110**

**For general laboratory use (RUO)**

**Store the kit at +15 to +25°C.**

**BENCH-TOP PROTOCOL**

---

Biotechnology Incubator, National Institute of Genetic Engineering and Biotechnology (NIGEB),  
Shahrak-e Pajoohesh, 15th Km, Tehran-Karaj Highway, Tehran, Iran

Contact us:

(+98) 21 44787435, (+98) 9033115921

[www.noyafan.ir](http://www.noyafan.ir)

## Purification of PCR Products

- 1** - Adjust total volume for each PCR tube to 100  $\mu$ l with water.  
- Add 500  $\mu$ l Binding Buffer per each 100  $\mu$ l PCR product. Mix sample (Binding Buffer + PCR solution) well.

---
- 2** - Place a DNA spin column in a provided 2 ml collection tube. Transfer the sample from step 1 to the upper reservoir of the DNA spin column.  
- Centrifuge 30 – 60 s at maximum speed at room temperature.

---
- 3** - Disconnect the filter tube, discard flow-through and place the filter tube back in the same tube.  
- Add 500  $\mu$ l Wash Buffer to the upper reservoir.  
- Centrifuge 1 min at maximum speed at room temperature. Discard the flow-through liquid.

---
- 4** - Add 200  $\mu$ l Wash Buffer and centrifuge as above. Discard the flow-through solution.  
Optional (recommended): Reconnect the filter tube to a clean 1.5 ml micro-centrifuge tube and centrifuge again to ensure full removal of Wash Buffer.

---
- 5** - Reconnect the Filter Tube to a clean 1.5 ml micro-centrifuge tube and add 50–100  $\mu$ l Elution Buffer directly to the center of the column membrane and stand for 3 min.  
- Centrifuge 1 min at maximum speed at room temperature.  
- Store eluted DNA at 4°C or -20°C.

---

Contact us:

(+98) 21 44787435, (+98) 9033115921

[www.noyafan.ir](http://www.noyafan.ir)

