

## Gel documentation system



The Paya- gel documentation system allows users to perform image capture functions with user-friendly panel. The simple design is excellent for multiple user laboratories. It can easily connect to computer by though USB cable to transfer images directly to your computer. The reliable and safe Paya gel documentation is an excellent choice for molecular biology laboratories.



### Technical specifications

<b>UV TRANSILLUMINATOR</b>	<b>256 nm &amp; 365nm</b>
<b>FILTER</b>	<b>580AF12062mm Ethidium Bromide Filter</b>
<b>Resolution</b>	<b>18 Megapixels</b>
<b>Lens</b>	<b>55-250mm</b>
<b>Data type</b>	<b>JPEG</b>
<b>UV Transilluminator size</b>	<b>20*20 cm</b>
<b>Electrical Requirement</b>	<b>220 V AC - 50HZ</b>
<b>WEIGHT</b>	<b>10 KG</b>
<b>DARKROOM DIMENSION (W*H*D)</b>	<b>34*67*44 CM</b>