



PCR enhancer

Cat No: DP9041

Store at: -20 °C

For research use only

In total 1 vial

Description

PCR Enhancer system is used in conjunction with DNA polymerases to optimize PCR of problematic and/or GC-rich templates. This system offers higher primer specificity, boarder magnesium concentration, boarder annealing temperature and improving themostability DNA polymerase.

Applications

Complex template with GC-rich

High efficiency, quality and yield

Quality control

Ad 2-8 µl PCR enhancer to 25 µl amplification system, the other amplification condition and procedure remain the same.

Product use limitation

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor suitable for administration to humans and animals.