


Inoculating Turntable Inoculation Turntable for Bacteria Culture Plates

IVD



Rotational speed 90 to 110 rpm
Accommodates Petri dishes up to 105 mm diameter On/Off sensor allows user to switch the turntable on and off at any time with simple movements of the hand.

A continuous rotational operation setting is for use when a longer operating period is required.
Dimensions W x H x D: 150 x 80 x 200 mm.



MST Park, Karmaco. <http://karmaco.epage.ir>

Fax: +9802144153205 Tehran-Iran

Karvandian@modares.ac.ir