

## Teksan Co.

Light Spectroscopy Instruments

### Address

Unit 104, Technology Units  
Incubator, Shahid Beheshti  
University, Tehran, Iran  
Postal code: 1983969411

### Contact information

Tel: +98 21 22402199  
Telegram: @Teksanco  
Fax: +98 21 43855749  
Linkedin.com/company/Teksan-co/

### Website and E-mail

www.teksan.ir  
Sales@teksan.ir  
info@teksan.ir

For more information on any of  
our products or services please  
visit us on the Web.

# Portable Raman (Portino - 532 nm Excitation)

Comprehensive, Fast and always with you

Full option and user-friendly control program



## Application Areas

- Materials Engineering
- Food & Beverage Quality Control
- Chemistry
- Pharmaceutic
- Biology
- Physics
- Protein & Nucleic Acid Analysis
- Volcanology
- Color Analysis
- Environments
- Nanotechnology
- Polymers
- Agriculture

The development of portable Raman spectrometers has drastically changed the way the industry conducts quality control for incoming raw materials. Raman has always provided advantages over other methods in that it does not require sample pre-treatment or direct contact with the sample, and has the unique capability of being able to test a sample directly through a transparent packing material like glass or plastic.

Portino is designed for rapid identification of materials, from field measurements of chemical and explosive agents to counterfeit detection and authentication analysis. Its ability to be used in harsh and demanding environments is well suited to quality assurance and verification of incoming materials, as well as in-line process monitoring.

- Fastest Sampling
- Large range of options
- All modes of microscopy, transmission and reflection illumination
- Extended Spectrum acquires
- Fiber optic coupling for convenient sampling
- Easy-to-use
- Portability for in-line, at-line, or laboratory analysis
- Customizable sample holders to accommodate any sample dimension
- Tight Optical Configurations
- Step by step, easy to follow quantitative modeling and method development