

# ET-Adsorb<sup>®</sup>

## Ecological Ethidium Bromide Removal System

The ET-Adsorb<sup>®</sup> Ethidium Removal System provides an elegant and easy way to remove large quantities of Ethidium from contaminated inorganic solutions. The cartridge has a total capacity of at least 1 g of Ethidium. This equals to 1,000 l of a typical electrophoresis buffer or staining solution (at a concentration of 0.5 µg/ml)

Saturated cartridges can be disposed of via incineration.



### Ideal for:

- ✓ Sustainable handling of Ethidium Bromide solutions
- ✓ High binding capacity: at least 1 g of Ethidium per cartridge, or 1,000 l of buffer or staining solution (at a concentration of 0.5 µg/ml)
- ✓ Easy to use
- ✓ Easy exchange of cartridges
- ✓ Buffer reservoir of 4 l
- ✓ Flow rate: approximately 2.0 l/hour
- ✓ Cartridge disposal via incineration

Kiagen Teb Co.

Unit31, 10<sup>th</sup> Floor, No 137, Dogol Tower, Opposite of Amir-Atabak  
Motahari St, Tehran

Tel: 021-88760179, 0912-3847794

Fax: 021-88759622

[www.kiagen-biotech.com](http://www.kiagen-biotech.com)

[info@kiagen-biotech.com](mailto:info@kiagen-biotech.com)