

# آزمایش اینتر فرومتر فابری پرو- مایکلسون

## interferometer Exp

### Required equipments for the Experiment:

- Sodium lamp with housing 110011
- Special stand 110002
- Power supply for sodium and mercury lamp 110010
- Fabry-Pro interferometer 220007
- High pressure mercury lamp with socket 110009
- Green filter for mercury lamp and holder 220004
- Stand base, V shape 110005
- Stand Rod 110007
- Saddle base 110006
- MultiClamp 110003
- Stand base V Shape (Small) 110005

### Descriptions:

Measuring the different wavelength between two spectral lines of Na (Yellow)  
Measuring the unknown wavelengths of light.

**Warranty :** One year

**After sales service :** Ten years

