

# neuroKit

نرو سیستم عرضه کننده انواع الکتروود و مواد مصرفی تحریک  
های الکتریکی

tDCS/HD-tDCS/tES



تماس با ما

تهران - خیابان باقرخان - بیمارستان امام خمینی - جنب مرکز پیوند اعضا -  
مرکز تحقیقاتی تجهیزات پزشکی پیشرفته - اتاق ۳۵

۰۲۱-۶۶۹۰۷۵۱۶



Neurosystem.ir




Electrode protection sponges are used to store saline material that can uniformly distribute the flow density throughout the electrode.

**ULTRA [ARRAY]**  
**PROCESSING SINA**

اسفنج های محافظ الکترود جهت ذخیره ماده سالین در خود بکار می روند که می توانند سبب توزیع یکنواخت چگالی جریان در سرتاسر الکترود گردند.



انواع اسفنج های تحریک با هر سایز و مطابق با نیاز مصرف کنندگان



## Electrode Size

PREMIUM

**Materials:** Conductive Rubber

**For 3x5 cm neuroPad**

**Dimensions:** 2.5 cm x 2.5 cm

**For 5x5 cm neuroPad**

**Dimensions:** 2.5 cm x 4.5 cm

**For 5x7 cm neuroPad**

**Dimensions:** 4.5 cm x 4.5 cm

**For 5x10 cm neuroPad**

**Dimensions:** 4.5 cm x 4.5 cm

neuroPad are designed to provide consistent performance over the duration of a tDCS/tES session, including setup time. Their hyper-absorbent material also allows supplementing saline during stimulation as needed, based on feedback from the SmartScan device sensor.

The only sponge electrode for tDCS/tES made from medical grade materials and certified bio-compatible (single use only). In medical settings, where safety and tolerability cannot be compromised, neurosystem neuroPad are used.

جهت سفارش الکترود خود می توانید از طریق مراجعه به سایت نروسیستم سفارش خود را ثبت نمایید.

