

# Teif Gostar Faraz

19 Inch 4 unit, WITS output

TGF

## Online Gas Chromatography

The OGC-2570 is the one of best choice of gas chromatography that analyses most gaseous component very fast. The major advantage for high quality analyses are the use of advanced GC configuration using multi column and multi port valves.

**Fast analysis of C1-C5 less than 100 second – THC report**

**Altitude compensation system**



**OGC-2570**

### Performance

High speed, High resolution, High reproducibility, High quality High sensitive FID and THC (Optional) or TCD

Double split capillary system split only C1 and C2 at mud logging mode

Accurate and high performance oven with any step temperature programming

Full EPC high resolution gas control for all inlet and detectors

**User friendly software and level zero WITS output.**


### Application

Permanent gas and Oil well logging gas analyses

Low weight gas analyses


Field and laboratory instrument

Natural gas calorific value calculation



Unit 4, No. 183, Mahtab Bldg.,  
Tolid Garan Blvd., KM 19  
Jade Ghadim Karaj,  
Tehran 3114859018, Iran

Tel: 0098 21 44336672  
Fax: 0098 21 46825835  
E mail: [info@irangc.com](mailto:info@irangc.com)  
Telegram: @irangc95  
Teif Gostar faraz Co. Ltd.



**طیف گستر فراز**

طراحی و ساخت گاز کروماتوگراف  
تهران: ۴۴۳۶۶۷۲