



# EEG 3520



24 Channel  
Electroencephalography

- ACCURATE
- RELIABLE
- EASY TO USE

● Full programmable functions such as:

- Montage manager
- Automatic pattern for EEG signals recording
- Automatic calibration

● Optional video EEG recording


● Strong trolley, and easy to move

● General specification:

- Low cut filter: programmable between 0.16 HZ- 3.2 Hz
- High cut filter: programmable between 10HZ- 90HZ
- Notch filter: 50 HZ
- Sample rate:250 Sample/Second
- Photic stimulator frequency: 0-20HZ

● The Montage Editor allows you to create new montage and modify previous montage.

● Selection program can help you to define a selected area and you can extract specific area from main EEG trace.



● The Navigator window can navigate you to the any selected area.

