NOYA DEPC-Treated Water

Catalog number: N-1181

«Research Use Only»

PRODUCT DESCRIPTION

18-megaohm ultrapure water treated with diethylpyrocarbonate (DEPC) and autoclaved to remove the residual DEPC. Sterile filtered and packaged in sterile bottles.

QUALITY CONTROL

RNAse testing: 7 μ g HCV RNA was incubated for 1 hour at 37°C with 10 μ l DEPC Treated water in a 100 μ l reaction consisting of 10mM NaCl, 5mM Tris-HCl pH 7.9, 1mM MgCl2, and 0.1mM dithiothreitol. No degradation of the RNA observed by real time PCR.

Size: 25mL

SHIPPING CONDITIONS: Ambient

STORAGE CONDITIONS: 2-8 ° C